



2.2mm

T

his blister is recommended to inject Micro lenses in acrylic hydrophilic hydrophobic

This kit is composed of :

- 1 cartridge MDJ Cart + reference 00265 Mini B (Round Bore Cartridge and System to clips winglets)
- 1 single used "Mini-Glider" injector with a silicone tip for a hydraulic injection
- 1 handling pusher to well position the IOL into the cartridge.

Use by date : 18 months

Label with 4 detachable stickers, slides of protection.

2.5mm

This blister is recommended to inject monobloc lenses in acrylic hydrophilic hydrophobic or silicone for in

This kit is composed of :

- 1 cartridge MDJ Cart + reference 00263 XSBI (Round Bore Cartridge and System to clips winglets)
- 1 single used "Mini-Glider" injector with a silicone tip for a hydraulic injection
- 1 handling pusher to well position the IOL into the cartridge.

Use by date : 18 months

Label with 4 detachable stickers, slides of protection.

2.8mm

T

This kit is composed of :

- 1 cartridge MDJ Cart + reference
- 1 single used "Mini Glider" injector with a silicone tip for a hydraulic injection
- 1 handling pusher to well position the IOL into the cartridge.

Use by date : 18 months



Label with 4 detachable stickers, slides of protection.

